



TECHNOLOGY

Microtentacle Imaging in Patient Tumor Samples

OVERVIEW

Free floating microenvironments are relevant to many aspects of healthcare and research. However using currently available technologies, it is difficult to utilize non-adherent cells for research and diagnostic purposes due to harsh processing steps (e.g. centrifugation, washing) or artificial substrates (e.g. plastic, glass) that change cell behavior.

This technology is a diagnostic tool to observe and study non-adherent cells in a free-floating environment that will preserve the biological behaviors of non-adherent cells. More specifically, this technology is a polyelectrolyte multilayer (PEM) surface coating designed with specialized lipids to tether cells to microfluidic substrates such as microscope slide platforms and cell culture surfaces.

The PEM coated slides have been used to image cell membrane protrusions (microtentacles) that form on the surface of free-floating tumor cells, allowing for the precision image capture of dynamic membrane structures. Microtentacles aid tumor cell reattachment to blood vessel walls and retention of circulating tumor cells (CTCs) in lung capillaries during metastasis. To date, this technology has been used to successfully conduct non-adherent cellbased analysis and screening on cancer cell lines and patient CTC samples.

APPLICATIONS

The primary application for this technology is as a live cell capture and imaging microfluidic device (i.e. microscopy platform slide or lab-on-a-chip) to examine the behavior of cells in a non-adherent, free-floating state. This would apply to native non-adherent cells (i.e. immunocytes) as well as rare cells or scarce samples such as CTCs, stem cells, and other anchorage-independent cancer cells.

ADVANTAGES

Microfluidic platform: Ability to immobilize individual or populations of live cells for analysis using a label-free, low-cost platform

-MEW

LICENSING POTENTIAL

Available for licensing in other fields other than oncology and research tools.

CONTACT INFO

Office of Technology Transfer
620 W Lexington St., 4th Floor
Baltimore, MD 21201
Email: ott@umaryland.edu
Phone: (410) 706-2380

Additional Information

INSTITUTION

University of Maryland, Baltimore

PATENT STATUS

SM-2013-121- EU 15 772 640.7 filed 11/2/2016, US 15/301,587 filed 10/3/2016 SM-2015-011 – PCT/US2015/054571 filed 10/7/2015

LICENSE STATUS

Available

INVESTIGATOR(S)

Stuart Martin

ATTACHMENTS

-  [Download SM-2013-121 SM-2015-011.pdf](#)

EXTERNAL RESOURCES

- [Lipid tethering of breast tumor cells enables real-time imaging of free-floating cell dynamics and drug response](#)

SM-2013-121. SM-2015-011