

TECHNOLOGY Membrane Displacement Traps for Droplet Discretization and Manipulation

OVERVIEW

Background

Microfluidics is a class of technology where the precise manipulation of fluids in very small channels and devices leads to improved devices in fields as different as medical diagnoses and inkjet printing. A subclass of microfluidic technology is Droplet-based Microfluidics where instead of manipulating a continuous flow, fluids are broken into discrete volumes (droplets) allowing diverse applications in chemical and biological analysis and processing such as the production of microcapsules and microparticles, digital PCR, protein crystallization, and single-cell analysis. Unfortunately, most droplet generating platforms are designed with complex architectures and operation mechanisms requiring ancillary equipment to achieve droplet discretization, metering, capture and coalescing. This becomes a hindrance for the commercialization of droplet-based microfluidic devices since it requires skilled operating technicians with little room to scale up.

Invention:

Researchers from the University of Maryland's Department of Mechanical Engineering have invented a novel droplet generation and manipulation platform for microfluidic devices that is simple to produce and easy to scale from single droplets to thousands of droplets. Featuring an array of individually addressable membrane displacement traps, this platform provides flexible and programmable control over droplet manipulation using a single actuator topology to perform all key droplet operations including discretization, release, transport, capture, metering, and merging. The platform is also designed to allow for bidirectional movement of droplets that will allow for arbitrary combinations of both serial processing steps.

Advantages:

- · Single actuation mechanism for multiple functions
- · Readily adapted for multiple analyses
- · Cost efficient due to reduction of operating equipment

Applications:

- · Nucleic Acid based assays
- · Single cell analyses
- · Drug Screening
- · Microparticle production and manipulation
- · Research tools

Research Article: Lab Chip, 2017, 17, 3717-3724 CONTACT INFO

UM Ventures 0134 Lee Building 7809 Regents Drive College Park, MD 20742 Email: <u>umdtechtransfer@umd.edu</u> Phone: (301) 405-3947 | Fax: (301) 314-9502

Additional Information

INSTITUTION

University of Maryland, College Park

EXTERNAL RESOURCES

PS-2017-121